

[Current Issue](#) [Previous Issues](#) [Contact](#) [Subscribe!](#)

Getting EEG data in real-time for BCI, Neurofeedback and more

Brain Products / [December 17, 2015](#)

by Dr. Mario Bartolo

Technical Consultant (Brain Products)

In the past years there has been an increased need and demand for real-time access to EEG data for purposes of quality control, parametrization, features extraction and online feedback. In this article, I will explore the concept of remote data access (RDA) and demonstrate how to access EEG data in real-time from our amplifiers.

› Overview

Not long ago EEG analysis were left running overnight in an attempt to overcome the long and cumbersome processing time that computers offered. Nowadays, EEG data analyses take a few milliseconds, so researchers started to “listen” to the human brain and “change” cognitive states/functions during the experiment.

Taking advantage of the high temporal resolution provided by EEG and the possibility to have online access to the data, BCI and Neurofeedback research fields rapidly developed and grew. **Brain Products always allowed researchers to have real-time access to data acquired with Brain Products amplifiers.** We offer this possibility on two levels:

- Software Development Kits (**SDK**) – access on driver level our amplifiers without using a recording software
- Remote Data Access (**RDA**) – stream EEG data to other programs while it is being recorded

The SDK option is out of the scope of this article and you can find detailed information on how to access the USB drivers (amplifier dependent) in our [website](#). **In this article, I will describe**

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

[Close](#)[More Information](#)

RDA is the term given to remotely access our recording software (**BrainVision Recorder**) and transfer data to other programs located on the local computer or on computers within the same network. In this process, the recording software acts as the server, and the program receiving the data acts as the client (e.g. BrainVision RecView). Data are transmitted via TCP/IP (Figure 1). It is important to highlight that it is possible to log up to ten clients to the server computer, whilst data are being saved at the server's hard drive.




Fig 1: Schematic representation of the interaction between server and client computer.

› Establishing Remote Data Access

In what follows, I will provide an example of how to establish a TCP/IP communication between a RDA Server Computer i.e. the Recorder computer and a RDA Client Computer where a Python script for online analysis is running. In this example the TCP/IP communication between two computers is established by using an Ethernet cable. It is important to highlight that if your computers are not equipped with Gigabit Ethernet, then you will need an Ethernet crossover cable.

Configuring the Server Computer




Fig. 2: RDA activation in BrainVision Recorder.

As stated above, the server computer refers to the computer where Recorder is running. Thus, the initial step is to enable the RDA option in Recorder. To do this, after opening Recorder, go to **Configuration › Preferences ...**, select the **Remote Data Access** tab and tick the **Enable Remote Data Access** box (Figure 2).

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

Close

More Information

Once the connection is established, please go to **Control Panel > Network and Sharing Center** and open the **Network Connection Details** of the newly established Ethernet Connection. Here, you will find the IPv4 address of the server computer (Figure 3) which we will use latter when configuring the client computer.

A last step is to enable communication to the server’s RDA service. More specifically, we need to ensure that the port used by Recorder’s RDA server (i.e. port 51244) is not blocked by any of Windows security features and/or installed security programs (e.g. anti-virus, firewall, Windows Defender ...). There are two possible ways to achieve this: either temporarily disabling all security programs or ensuring that the port used by Recorder’s RDA server is added to the exceptions list of your security program(s). While the former option is simpler and convenient as a temporary solution, the latter is more complex to achieve but nevertheless safer as a long term solution.




Fig. 3: Network Connection Details: Server computer IP address highlighted in blue.

Configuring the Client Computer

The first task of the client program is to establish a connection to the server’s RDA service. This is done using standard socket programming which will not be covered here. Instead, I will describe how to configure one of our RDA clients’ scripts that are freely available in our [website](#): RDA Client for Python (RDA.py).

Figure 4 shows an excerpt of the above mentioned Python script. Here, the red rectangle is highlighting the code line where the server address and port number needs to be defined. From the previous section (Figure 3) the server address (where Recorder is running) was identified as being 169.254.7.107. Please be aware that the server address might differ every time a new connection is established and so the code line needs to be modified accordingly.



Fig. 4: Excerpt of the RDA Client for Python (RDA.py). The red rectangle identifies the code line where to add the server IP address and port.

RDA Client for Python output

In the last two sections, it was described how to configure the server and the client computer. The next step would be to start Recorder in Monitoring mode and run the RDA.py at the client computer (please ensure that you have Python software installed at the client computer). Once the connection is established, the client then waits for data or messages to be sent from the server. Further details about the contents and structure of the transmitted data block can be found in the [Recorder manual](#), more specifically in the Remote Data Access chapter.

Figure 5 shows the console output at the client computer when running the RDA Client for Python. In this specific example, and among various output messages, the average power of all channels is calculated over a time period of one second and printed into the console.

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

[Close](#)[More Information](#)




Fig. 5: Windows command line output when running the RDA Client for Python.

› Drawbacks and Limitations of RDA

TCP/IP communication is well known to be a reliable block-based transmission protocol. Thus, we have adopted it in order to allow researchers to remotely access their data in real-time. But any decision comes with a cost, and we are aware that when accessing the data remotely, a delay in the real-time transmission should be expected.

Taking into account the overall EEG signal pipeline (i.e. amplifier hardware filters → A/D converter → USB communication → Windows OS → RDA server (Recorder) → Ethernet Connection (TCP/IP) → RDA Client), average delay times for RDA are in the 50 – 100 ms range. This delay strongly depends on the RDA signal pipeline, heavy processing time, and usage of TCP/IP connection by another application. Even though in most applications the delay time is within the acceptable range, it is indeed one of the RDA limitations.

That said, RDA offers the most convenient, but not the fastest route to access data in real-time. Thus, and if you work with real-time critical applications, please keep in mind that there is always the possibility to access directly the USB drivers of our amplifiers without using a recording software (please see above. SDK) and therefore greatly decrease the delay time.

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

Close

More Information

Accessing data in real-time from Brain Products amplifiers is possible either by RDA or SDK. The existing **RDA option gives you full flexibility to program your online application in your favourite and preferred scripting language**. Furthermore, and in order to ease the process of implementing and programming your online analysis programs, Brain Products offers free access to **sample codes of RDA clients** (unsupported) in the most common programming languages (Matlab, C#, Python, C++). These sample codes could be used as a starting point/programming template for your own application.

©Brain Products GmbH 2015

...

December 17, 2015 in 2015, Issue 4/2015, Support & Tips. Tags: BCI, Brain Computer Interfaces, BrainVision Recorder, Neurofeedback, RDA, real-time, Remote Data Access, Support & Tips

Related posts

News in brief

BrainVision Analyzer 2 Webinars ... because you are never too far away to attend an Analyzer training session

Recording brain waves in extreme situations: the psychoactivity of ayahuasca

[Release of BrainVision Analyzer 2.1.1; usability, another step along the path](#)

[Brain Products Conferences and Events – January to May 2016](#)

Top Articles

All you ever wanted to know about markers in BrainVision Analyzer 2

Independent Component Analysis - demystified!

Track the eye while EEG is acquired: Now also with SR EyeLink 1000 Plus Eye Tracker

Getting EEG data in real-time for BCI, Neurofeedback and more

Articles by category

Select Category



Search the “Press Release”

All Categories



Search

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

Close

More Information

Tweets by @Brain_Products

Brain Products Retweeted

**Emily Kappenman**

@EmilyKappenman


A full-time research assistant position is available in my cognitive/affective neuroscience lab starting this summer! For more information, contact me at emily.kappenman@sdsu.edu
Please share/retweet

Apr 26, 2019

**Brain Products**

@Brain_Products

Visit our partner, INTER, at the 6th Intl. Symposium on Brain and Cognitive Science (ISBCS) on April 28th in [#Istanbul](#), Turkey. Stop by the INTER booth to learn more about Brain Products solutions for [#EEG](#) [#research](#). [@YUCognitive](#) sta.cr/36xi1



SIAIS IS THE SATELLITE EVENT
OF INTERNATIONAL SYMPOSIUM
ON BRAIN AND COGNITIVE

Apr 26, 2019

**Brain Products**

@Brain_Products

We just published another [#Analyzer](#) 2 [#webinar](#) recording. Visit the Brain Products Webinar Channel and watch "Introduction to Source Analysis - Distributed Source Imaging using [#LORETA](#)" at any time. sta.cr/36Ed1



This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

Close

More Information

Apr 25, 2019

[Embed](#)

[View on Twitter](#)

Proudly powered by [WordPress](#) | Theme: Expound by [Konstantin Kovshenin](#)

5

This website uses cookies to improve your experience. By using and further navigating this website you accept this. Detailed information about the use of cookies on this website is available by clicking on more information.

Close

More Information